

# RESCUE 230

*Defibrillatore Bifasico*

*Biphasic Defibrillator*



**Immediate  
Action  
and simple  
to use**



- BTE (biphasic truncated exponential) discharge type
- Impedance compensated energy delivery
- Energy selection from 10 to 230 Joule directly from the front panel.
- Charge time 5 sec for 230 Joule with fully charged battery
- 100 shocks of 230 J with fully charged batteries
- Charge/discharge control from the paddles
- Ecg 'R' wave synchro mode or async mode
- Can be used with disposable paddles
- Large (5.7") high resolution backlite LCD display (320x240 dots)
- ECG through the paddles or through patient cable

***An Italian Product***

- BTE (onda bifasica esponenziale troncata)
- Scarica dell'energia controllata tramite misurazione di impedenza
- Selezione di energia da 10 a 230 Joule direttamente dal pannello
- Tempo di carica: 5 secondi per 230 Joule con le batterie cariche
- 100 scariche da 230 J con batterie completamente cariche
- Carica/Scarica controllabili direttamente dai manipoli
- Modo sincrono ed asincrono
- Può essere usato con piastre monouso
- Display LCD (5.7") retroilluminato alta risoluzione (320x240 punti)
- Rilevamento ECG tramite i manipoli di scarica o il cavo paziente

***Rapidita' e Semplicita' d'uso***



**PROGETTI**

MEDICAL EQUIPMENT